

CF-LAB/Active Carbon. Central Extractor Filter.

Active Carbon 30 Ltr.

CF-LAB/Active Carbon Central Extractor Filter is equipped with a 30Ltr Active Carbon Filter Cartridge with a Pre-Filter Fleece on the outside. It is suitable for filtering gases and odours.

Application Examples.

1. Cleaning / Degreasing
2. Conservation
3. Gluing / Bonding
4. Moulding
5. Screen Printing
6. Mixing, Measuring, Taking Samples, Examinations etc.

Replacement Filter Cartridge.

It is time to put in a Replacement Filter when the extraction is not satisfactory or the odours come back. It is easy to change the Filter Cartridge and will just take a few minutes.

Capacity.

4-8 Flexa/Chemikus Extraction Arms (d. 75-100mm) can be connected to a CF-LAB/Active Carbon Central Extractor Filter or 1-2 Flexi Extraction-Arms (d.160mm/max.2m). The inlet diameter is 199mm. Use d.200mm Vertical/Horizontal ducting from the CF and go down with 75/100mm ducts to the Flexa/Chemikus Extraction Arms and d.160mm duct to Flexi dia.160mm Extraction-Arm/Arms.



Ask our specialists for further advice.

40 Years Experience!

40 Years Experience of design and manufacturing at your service!

The Plymoth products are a combination of well-proven technique and progressive new ideas.

Plymoth is always a step ahead - but not expensive!

	Art. No.
CF-LAB/Active Carbon. ~1/1.1kW/220V/50Hz. 30Ltr. Act.Carb. Inlet 199mm. Legs (0.5m).	P-1012
CF-LAB/Active Carbon. ~3/1.1kW/400V/50Hz. 30Ltr. Act.Carb. Inlet 199mm. Legs (0.5m)	P-1013
Replacement Active Carbon Cartridge 30Ltr. For MF, WF, CF, MF-Filter Tables (High).	P-962

J. Plymoth AB, Trarydsallén 4E, 287 72 TRARYD, Sweden.

Tel.: +46 40 700 15 · E-Mail: info@plymoth.com · www.plymoth.com